

Date: 03/02/2021

Using the NObreath[®] during COVID-19 (1st and 2nd generation)

Key organizations such as the Primary Care Respiratory Society (PCRS) and the Association of Respiratory Technology and Physiology (ARTP) have advised that FeNO testing is **low risk** and **unlikely to generate infectious aerosols**, therefore supporting FeNO testing as a '**Non-Aerosol Generating Procedure**' (Non-AGP).

It is also noteworthy that there have been no reports of coughing experienced from the use of any of Bedfont[®] Scientific's first and second generation NObreath[®] device. Recent clinical and usability studies performed on our NObreath[®] devices reviewing safety and performance, **have no reported adverse events** where coughing was seemingly induced by the breathing maneuvers required by Bedfont[®] Scientific's NObreath[®] device.

Despite the safety of FeNO testing in general; both the 1st and 2nd generation NObreath[®] mouthpieces contain an infection control filter, which has been independently tested for efficacy by Public Health England. The filtration rates are displayed below for both 1st and 2nd Generation NObreath[®] mouthpieces:

1 st Generation NObreath [®] Mouthpiece Bacterial Filtration Efficiency (BFE) and Viral Filtration Efficiency (VFE)	
BFE % rate	>99%
VFE % rate	>99%

2 nd Generation NObreath [®] Mouthpiece Bacterial Filtration Efficiency (BFE) and Viral Filtration Efficiency (VFE)	
BFE % rate	>99%
VFE % rate	>98%

Both 1st and 2nd generation NObreath[®] mouthpieces have been tested to filter viruses as small as 24 nanometres in diameter and the Covid-19 virus particle has a diameter of approximately 125 nanometres. The virus model used is incredibly penetrable, even more so than a majority of human viruses, therefore makes it a very effective model to use for virus filtration efficiency (VFE). Therefore, Bedfont[®] can conclude that bacterial and viral pathogens (including COVID-19) will effectively be removed by the 1st and 2nd Generation NObreath[®] mouthpiece filter at the efficiency rates shown above.

The NObreath[®] mouthpieces are a single-patient use mouthpiece, meaning it should be disposed of according to local waste guidelines immediately after testing is complete to further minimize the risk of cross infection.

Furthermore, When taking a FeNO measurement with the NObreath[®], due to our uniquely designed NO scrubber and software algorithms, the patient does not inhale through the device or mouthpiece before exhaling, to further reduce the risk of cross-infection.



Jason Smith
Managing Director

Bedfont Scientific Ltd

Station Road, Harrietsham, Maidstone, Kent, ME17 1JA, England
Tel: +44 (0)1622 851122, Fax: +44 (0)1622 854860, Email: ask@bedfont.com
Registered in: England and Wales. Registered No: 1289798 Issue 3 - February 2017, Part No: MKT361