

NObreath®

Mouthpiece Datasheet



Designed specifically for the NObreath® FeNO monitor, the NObreath® mouthpiece contains an integrated infection control filter that removes and traps >99% of airborne bacteria and >98% of viruses¹.

Technical specification

Dimensions	Approx. 180 x 28 x 22 mm
Weight	Approx. 14g
Material	Polypropylene



The NObreath® mouthpiece may be used multiple* times by one patient, per clinical visit.

NOTE: The mouthpiece is single-patient use; re-use of the mouthpiece between patients can increase the risk of cross-infection and result in incorrect readings.

The NObreath® mouthpiece is fully recyclable and should be disposed of after use in accordance with local clinical waste guidelines.

References:

1. Public Health England. An Evaluation of Filtration Efficiencies Against Bacterial and Viral Aerosol Challenges. Salisbury: Public Health England; 2020.

*Single patient use for up to 15 tests per clinical visit.



Bedfont® Scientific Ltd.
Station Road, Harrietsham, Maidstone,
Kent, ME17 1JA England
Tel: +44 (0)1622 851122 Fax: +44 (0)1622 854860
Email: ask@bedfont.com Web: www.bedfont.com



Bedfont® GmbH
Mitterbachweg 18, 5081, Anif, Austria
Tel: +43 662 635050 Fax: +43 662 635050-50
Email: bedfont-gmbh@bedfont.com

Issue 3 - December 2020, Part No: MKT514

Bedfont® Scientific Limited reserves the right to change or update this literature without prior notice.

Registered in: England and Wales. Registered No: 1289798

